AM-SDG Medium 45 %
*ARTSMedia Sperm Density Gradient Medium 45 %***AM-SDG Medium 90 %**
*ARTSMedia Sperm Density Gradient Medium 90 %**HEPES buffered medium*
Medio tamponado con HEPES**Product No.**
N.º de productoAM-SDG Medium 45 %:
503020, 503050, 503100AM-SDG Medium 90 %:
504020, 504050, 504100**Technical Support**
Soporte técnico support@artsmedia.dk +45 5350 4530 www.artsmediadenmark.comARTSMedia Denmark
Kongevejen 149,
DK-2830 Virum, Denmark**Indications for use**

AM-SDG Medium 45% and AM-SDG Medium 90% is a ready to use two layer gradient system designed for separating human sperm using density gradient centrifugation.

The AM-SDG media system is suitable for semen preparation for Intra Uterine Insemination (IUI), In Vitro Fertilization (IVF), and IntraCytoplasmic Sperm Injection (ICSI).

Device Description

20 ml, 50 ml, and 100 ml glass bottles, individually packed, provided sterile. The AM-SDG media system is a HEPES-buffered medium with suspended silane-coated colloidal silica particles.

For professional use only.

Quality Control Specifications

pH: 7.20 – 7.90
(Release criteria: 7.20-7.60)

Osmolality:
AM-SDG Medium 45 %:
310-340 mOsm/kg
AM-SDG Medium 90 %:
320-350 mOsm/kg

Density:
AM-SDG Medium 45 %:
1.0500-1.0700 g/ml
AM-SDG Medium 90 %:
1.1050-1.1150 g/ml

Sterility: No growth
(Ph Eur. 2.6.1, USP <71>)

Endotoxins: < 0.5 EU/ml (USP <85>)
Human sperm survival assay
(% motility compared with control after 4 hours): ≥ 80%
Human sperm survival assay
(% motility compared with control after 24 hours): ≥ 75%

Use of Ph Eur or USP grade products if applicable

A Certificate of Analysis, for each batch is available at
www.artsmediadenmark.com

Storage

Store in original container between 2-25°C and avoid exposure to sunlight. Do not freeze. The products can be used within 7 days of opening when maintained under sterile conditions and stored at 2-8°C. The products remain stable after transport (up to a maximum of 5 days) at elevated temperatures (≤ 37 °C).

Warnings

- Do not use after expiry date.
- Do not use if the product is cloudy or shows any evidence of microbial contamination.
- Do not use if seal of the container is damaged or open upon delivery.
- Do not re-sterilize after opening.
- Discard excess (unused) media that has been warmed in accordance with local regulations for disposal of medical devices.

Precautions

- Caution: U.S. Federal law restricts this device to sale by or on the order of a physician.
- Media use should be conducted using aseptic techniques.
- Handle all samples as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling samples.

Instruction for use

AM-SDG Medium 45 % and 90 % should be used according to the established procedures and protocols optimized for the individual clinic.

Bring all media and centrifuge tubes to room temperature before use.

Semen preparation

- 1 Transfer 2.5ml AM-SDG Medium 45 % into the centrifuge tube.
- 2 Gently place 2.5ml AM-SDG Medium 90 % under the upper layer. Take care that the two layers are distinctly separated. The two-layer gradient is stable for up to two hours.
- 3 Gently place the liquefied/frozen semen sample onto the upper layer.
- 4 Centrifuge for 18 minutes at 350 g to 400 g.
- 5 Remove supernatant.
- 6 Add 2-3ml Semen Wash medium and resuspend the pellet.
- 7 Centrifuge for 10 minutes at 300 g.
- 8 Remove supernatant and repeat steps 6 and 7.
- 9 Remove supernatant and resuspend in Semen Wash medium.

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