AM Oocyte Aspiration

ARTSMedia Oocyte Aspiration Medium

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ARTSMedia Denmark
Kongevejen 149,
DK-2830 Virum, Denmark

Indications for use

AM Oocyte Aspiration is intended for retrieval of human oocytes by follicular aspiration or flushing. The medium can also be used for washing of oocytes.

Device Description

20 ml, 50 ml and 100 ml glass bottles, individually packed, provided sterile. AM Oocyte Aspiration is a HEPES-buffered balanced salt solution, with Physiological Salts Gentamicin Sulphate 10 µg/ml, Sodium Pyruvate, Sodium Bicarbonate, and Human Serum Albumin (HSA).

For professional use only.

Quality Control Specifications

pH: 7.2-7.6 (Ph Eur. 2.2.3, USP <791>)
Osmolality: 270-300 mOsm/kg
(Ph Eur. 2.2.35, USP <785>)

Appearance: Clear and particulate free
Sterility: No growth
(Ph Eur. 2.6.1, USP <71>)

Endotoxins: < 0.25 EU/ml (USP <85>)
1-cell Mouse Embryo Assay (MEA):
≥ 80% of embryos developed to expanded blastocysts at 96 hours following 60 minutes exposure to ARTSMedia Oocyte Aspiration Medium.

A certificate of Analysis, for each batch is available at www.artsmediadenmark.com

Storage

Store in original container between 2-8 °C and avoid exposure to sunlight. Do not freeze.

Warnings

- Do not use after expiry date.
- Do not use if the product is cloudy or shows any evidence of microbial contamination.
- Do not use if seal of the container is damaged or open upon delivery.
- Do not use medium that shows any evidence of particulate matter.
- Do not re-sterilize after opening.
- Do not re-use.
- Single-use only. Discard excess (unused) media that has been warmed in accordance with local regulations for disposal of medical devices.

- AM Oocyte Aspiration contains gentamicin sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic (allergic reaction).

Precautions

- Caution: U.S. Federal law restricts this device to sale by or on the order of a physician.
- Media use should be conducted using aseptic techniques.
- Always wear protective clothing when handling samples.
- This product contains albumin, a derivative of human blood. All blood products should be treated as potentially infectious. Source material used for this product was tested and found negative for HBsAg, Anti-HIV-1/-HIV2, HCV, HCV-RNA, and HIV-RNA. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated.

General notes before use

AM Oocyte Aspiration is HEPES buffered and should be pre-heated to 37 °C prior use.

Instruction for use

AM Oocyte Aspiration should be used according to the established procedures and protocols optimized for the individual clinic.

Oocyte collection

During oocyte collection, maintain AM Oocyte Aspiration at 37 °C.

Oocytes are highly sensitive to pH alterations and temporary cooling. Even slight temperature fluctuations can lead to disruption of the meiotic spindle, potentially causing chromosome dispersal. It's important to ensure the cap of the bottle is securely closed during the pre-warming process. Use a rubber-free syringe to draw approximately 20-50 ml of the solution at a time. Follow the manufacturer's instructions for assembling the needle and tubing set.

Fill the syringes with AM Oocyte Aspiration for flushing the follicle. Aspirate the follicular fluid and if flushing, flush the follicle using the same volume of medium as the amount of follicular fluid collected. Be mindful not to enlarge the follicle beyond its original size.

Carefully wash the oocytes using AM Oocyte Aspiration or pre-equilibrated culture medium, ensuring the removal of any blood clots and granulosa cells.

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