## AM Semen Wash

### ARTSMedia Semen Wash Medium

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ARTSMedia Denmark  
Kongevejen 149,  
DK-2830 Virum, Denmark

### Indications for use

AM Semen Wash is intended for washing of sperm and for sperm swim-up procedures. The medium can also be used for IUI (Intra Uterine Insemination).

### Device Description

20 ml, 50 ml and 100 ml glass bottles, individually packed, provided sterile. AM Semen Wash is a HEPES-buffered balanced salt solution, with Physiological Salts, Gentamicin Sulphate 10 µg/ml, Sodium Pyruvate, Sodium Bicarbonate and Human Serum Albumin (HSA).

*For professional use only.*

### Quality Control Specifications

pH: 7.2-7.6  
(Ph Eur. 2.2.3, USP <791>)  
Osmolality: 270-300 mOsm/kg  
(Ph Eur. 2.2.35, USP <785>)  
Appearance: Clear and particulate free  
Sterility: No growth  
(Ph Eur. 2.6.1, USP <71>)  
Endotoxins: < 0.25 EU/ml (USP <85>)  
Human sperm survival assay (HSSA)  
(% motility compared with control after 24 hours): ≥ 80%

A Certificate of Analysis, for each batch is available at [www.artsmediadenmark.com](http://www.artsmediadenmark.com)

### Storage

Store in original container between 2-8 °C and avoid exposure to sunlight. Do not freeze.

### Warnings

- Do not use after expiry date.
- Do not use if the product is cloudy or shows any evidence of microbial contamination.
- Do not use if seal of the container is damaged or open upon delivery.
- Do not use medium that shows any evidence of particulate matter.
- Do not re-sterilize after opening.
- Do not re-use.
- Single-use only. Discard excess (unused) media that has been warmed in accordance with local regulations for disposal of medical devices.

- AM Semen Wash contains gentamicin sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic (allergic reaction).

### Precautions

- Caution: U.S. Federal law restricts this device to sale by or on the order of a physician.
- Media use should be conducted using aseptic techniques.
- Always wear protective clothing when handling samples.
- This product contains albumin, a derivative of human blood. All blood products should be treated as potentially infectious. Source material used for this product was tested and found negative for HBsAg, Anti-HIV-1/-HIV2, HCV, HCV-RNA, and HIV-RNA. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated.

### General notes before use

AM Semen Wash is HEPES buffered and should be pre-heated to 37 °C prior use.

### Instructions for use

AM Semen Wash should be used according to the established procedures and protocols optimized for the individual clinic.

When performing IUI, follow the instructions of the specific catheters used.

### Washing of spermatozoa - example:

Spermatozoa can be washed at either room temperature or at 37 °C.

- Add 5 ml of AM Semen Wash to the liquified semen sample and completely mix the solution.
- Centrifuge the mixture for 5-10 minutes at approximately 300g.
- Carefully remove the supernatant, ensuring that approximately 0.5 ml of semen remains in the centrifuge tube.
- Add 5 ml of AM Semen Wash to the test tube and gently mix the solution until the pellet is fully dissolved.
- Centrifuge the solution again for 5-10 minutes at 300g.
- Remove the supernatant and resuspend the pellet in an appropriate volume of AM Semen Wash.

- The semen sample is now ready for IUI or for further use in the IVF process.

### Sperm swim-up procedure - example:

- Gently stratify 1 ml of AM Semen Wash above the prepared semen sample in a conical tube.
- Incline the tube at an angle of 45° and incubate for 30 to 60 min at 37 °C.
- Carefully collect the supernatant, which contains the motile sperm and transfer to a clean tube.
- Dilute with 5 ml AM Semen Wash and centrifuge the mixture for 5-10 minutes at approximately 300g.
- Remove the supernatant and resuspend the pellet in an appropriate volume of AM Semen Wash
- The semen sample is now ready for IUI or for further use in the IVF process.

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