**AM EmbryoGrip**

ARTSMedia Hyaluronan Embryo

Transfer Medium

Bicarbonate buffered

Bicarbonato tamponado

**Product No.**

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511010

510010

**Technical Support**

Soporte técnico

✉ support@artsmedia.dk

☎ +45 5350 4530

🌐 www.artsmediadenmark.com

ARTSMedia Denmark

Kongevejen 149,

DK-2830 Virum, Denmark

**Indications for use**

AM EmbryoGrip is intended for embryo transfer.

**Device Description**

10 ml glass bottles, individually packed, provided sterile. AM EmbryoGrip is a bicarbonate buffered balanced salt solution, with high concentration of hyaluronan, minerals, vitamins, amino acids, synthetic serum replacement, human serum albumin (HSA), gentamicin sulphate 10 µg/ml, and phenol red (only 510010).

**For professional use only.****Quality Control Specifications**

pH: 7.2–7.45 (Ph Eur. 2.2.3, USP <791>)  
Osmolality: 257-275 mOsm/kg  
(Ph Eur. 2.2.35, USP <785>)

Appearance: Clear and particulate free  
Sterility: No growth  
(Ph Eur. 2.6.1, USP <71>)

Endotoxins: < 0.25 EU/ml (USP <85>)  
1-cell Mouse Embryo Assay (MEA): ≥ 80% of embryos developed to expanded blastocysts at 96 hours

A certificate of analysis is available at [www.artsmediadenmark.com](http://www.artsmediadenmark.com)

**Storage**

Store in the original container between 2-8 °C and avoid exposure to sunlight. Do not freeze.

**Warnings**

- Do not use after expiry date.
- Do not use if the product is cloudy or shows any evidence of microbial contamination.
- Do not use if the seal of the container is damaged or open upon delivery.
- Do not use medium that shows any evidence of particulate matter.
- Do not re-sterilize after opening.
- Do not re-use.
- Single use only.
- Discard excess (unused) media that has been warmed in accordance with local regulations for disposal of medical devices.
- AM EmbryoGrip contains gentamicin sulphate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic (allergic reaction).

**Precautions**

- **Caution:** U.S. Federal law restricts this device to sale by or on the order of a physician.
- Media use should be conducted using aseptic techniques.
- Always wear protective clothing when handling samples.
- This product contains albumin, a derivative of human blood. All blood products should be treated as potentially infectious. Source material used for this product was tested and found negative for HBsAg, Anti-HIV-1/-HIV2, HCV, HCV-RNA, and HIV-RNA. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated.

**General Notes Before Use**

AM EmbryoGrip is bicarbonate buffered and should be pre-equilibrated in a CO<sub>2</sub> incubator set at 37 °C and 5-6% CO<sub>2</sub>.

**Note:** Adjust incubator settings to a higher CO<sub>2</sub> concentration if the altitude is above sea level.

**Instructions for Use**

AM EmbryoGrip should be used according to the established procedures and protocols optimized for the individual clinic.

**Example**

Transfer the embryos to the AM EmbryoGrip medium, which should be equilibrated with CO<sub>2</sub>. This medium resembles the IVC and allows for longer storage of embryos in the incubator prior to transfer. It also contains a high concentration of hyaluronan to enhance embryo implantation. Since it is only bicarbonate buffered, avoid letting the embryos sit outside the incubator for more than a few moments.

Prepare the embryos for uterine transfer in 20-35 µl of AM EmbryoGrip.

Before use, flush the transfer catheter thoroughly with AM EmbryoGrip and transfer immediately after loading embryos to catheter to maintain optimal pH.

Inspect the catheter to confirm that all embryos have been successfully transferred.

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