AM Handling Medium

ARTSMedia Gamete Handling Medium

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ARTSMedia Denmark
Kongevejen 149,
DK-2830 Virum, Denmark

Indications for use

AM Handling Medium is intended for handling of oocytes and embryos outside the incubator. The medium can also be used for embryo transfer and biopsy.

Device Description

100 ml glass bottles, individually packed, provided sterile. AM Handling Medium is a bicarbonate and HEPES buffered balanced salt solution, gentamicin sulphate 10 µg/ml, sodium pyruvate and human serum albumin (HSA).

For professional use only.

Quality Control Specifications

pH: 7.2–7.6
(Ph Eur. 2.2.3, USP <791>)
Osmolality: 270-300 mOsm/kg
(Ph Eur. 2.2.35, USP<785>)
Appearance: Clear and particulate free
Sterility: No growth
(Ph Eur. 2.6.1, USP<71>)
Endotoxins: < 0.50 EU/ml (USP <85>)
1-cell Mouse Embryo Assay (MEA): ≥ 80% of embryos developed to expanded blastocysts at 96 hours (following 60 minutes exposure to AM Handling Medium).

A certificate of analysis is available at www.artsmmediadenmark.com

Storage

Store in the original container between 2-8 °C and avoid exposure to sunlight. Do not freeze.

Warnings

- Do not use after expiry date.
- Do not use if the product is cloudy or shows any evidence of microbial contamination.
- Do not use if the seal of the container is damaged or open upon delivery.
- Do not use medium that shows any evidence of particulate matter.
- Do not re-sterilize after opening.
- Do not re-use.
- Single use only.
- Discard excess (unused) media that has been warmed in accordance with local regulations for disposal of medical devices.
- AM Handling Medium contains gentamicin sulphate. Appropriate precautions should be taken to ensure that

the patient is not sensitized to this antibiotic (allergic reaction).

Precautions

- Caution:** U.S. Federal law restricts this device to sale by or on the order of a physician.
- Media use should be conducted using aseptic techniques.
- Always wear protective clothing when handling samples.
- This product contains albumin, a derivative of human blood. All blood products should be treated as potentially infectious. Source material used for this product was tested and found negative for HBsAg, Anti-HIV-1/HIV2, HCV, HCV-RNA, and HIV-RNA. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated.

General Notes Before Use

AM Handling Medium is HEPES buffered and does not require CO₂ equilibration.

Instructions for Use

AM Handling Medium should be used according to the established procedures and protocols optimized for the individual clinic.

Example

Transfer: Transfer the embryos to the well with AM Handling Medium preheated to 37 °C.

The embryos are then prepared and transferred to the uterus in 20-35 µl AM Handling Medium.

Flush the transfer catheter well with AM Handling Medium prior to use.

Examine the catheter to ensure all embryos have been transferred.

Denudation

Remove oocytes from hyaluronidase drops after no longer than one minute while using a denudation pipette.

Transfer the oocytes to a fresh dish with AM Handling Medium to wash away any excess hyaluronidase and dissociated cells. Gently pipette the oocytes in the fresh medium several times to ensure thorough removal of all residues, should any cumulus cells remain.

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