AM HyaSpermSelect

ARTSMedia Hyaluronan Sperm Select Medium

HEPES and bicarbonate buffered

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ARTSMedia Denmark
Kongevejen 149,
DK-2830 Virum, Denmark

Indications for use

AM HyaSpermSelect is intended for decreasing sperm movement, for the selection of the most mature and viable spermatozoa for ICSI.

Device Description

0.5 ml in plastic bottles, individually packed, provided sterile. AM HyaSpermSelect is a bicarbonate and HEPES buffered balanced salt solution, with high concentration of hyaluronan, minerals, synthetic serum replacement, human serum albumin (HSA), and gentamicin sulphate 10 µg/ml.

For professional use only.

Quality Control Specifications

pH: 7.2–7.6
(Ph Eur. 2.2.3, USP <791>)
Osmolality: 270-300 mOsm/kg
(Ph Eur. 2.2.35, USP <785>)
Appearance: Clear and particulate free
Sterility: No growth
(Ph Eur. 2.6.1, USP<71>)
Endotoxins: < 0.50 EU/ml (USP <85>)
Human sperm survival assay (HSSA) (% motility compared with control after 30 minutes): ≥ 80%

A certificate of analysis is available at www.artsmmediadenmark.com

Storage

Store in the original container between 2-8 °C and avoid exposure to sunlight. Do not freeze.

Warnings

- Do not use after expiry date.
- Do not use if the product is cloudy or shows any evidence of microbial contamination.
- Do not use if the seal of the container is damaged or open upon delivery.
- Do not use medium that shows any evidence of particulate matter.
- Do not re-sterilize after opening.
- Do not re-use.
- Single use only.
- Discard excess (unused) media that has been warmed in accordance with local regulations for disposal of medical devices.
- AM HyaSpermSelect contains gentamicin sulphate. Appropriate precau-

tions should be taken to ensure that the patient is not sensitized to this antibiotic (allergic reaction).

Precautions

- Caution:** U.S. Federal law restricts this device to sale by or on the order of a physician.
- Media use should be conducted using aseptic techniques.
- Always wear protective clothing when handling samples.
- This product contains albumin, a derivative of human blood. All blood products should be treated as potentially infectious. Source material used for this product was tested and found negative for HBsAg, Anti-HIV-1/-HIV2, HCV, HCV-RNA, and HIV-RNA. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated.

General Notes Before Use

AM HyaSpermSelect is bicarbonate and HEPES buffered and does not require CO₂ equilibration. It can be pre-equilibrated in a CO₂ incubator set at 37 °C and 5-6% CO₂. To maintain pH, temperature, and osmolality, make sure that the medium is completely covered with oil.

Instructions for Use

AM HyaSpermSelect should be used according to the established procedures and protocols optimized for the individual clinic.

Example

- Preheat oil and AM Handling Medium or if AM-IVC Medium is used for ICSI, pre-equilibrate in CO₂. The final prepared dish with oil overlay can also be equilibrated.
- Add 5-10 µL droplets for ICSI medium depending on the number of oocytes for injection.
- Add a drop of AM HyaSpermSelect, the volume depending on the laboratory's protocol. 5-10 µL for a small drop or elongated drop up to 200 µL.
- Cover the dish immediately with an oil overlay to prevent evaporation of media drops and incubate or keep at 37 °C degrees.
- Add one or multiple oocytes to the ICSI

drops according to laboratory protocol.

6. Add a small volume, 1-3 µL of sperm solution into the periphery of the AM HyaSpermSelect droplet.

7. Select one or several spermatozoa for injection and break the tails of the spermatozoon with the tip of the pipette. Spermatozoa with decreased motility indicate they are mature binding to the hyaluronan.

8. The spermatozoa can now be transferred through the oil to the ICSI drops and oocytes injected.

9. Transfer the injected oocytes to culture dishes with CO₂ equilibrated AM-IVC Medium for culture.

Note: Mature sperm possess hyaluronan receptors, allowing them to bind to hyaluronan in this medium. This binding results in increased viscosity, causing the sperm to swim slower. Therefore, selecting these sperm for Intracytoplasmic Sperm Injection (ICSI) can enhance the chances of successful fertilization by ensuring that only the most mature, functional sperm are used.

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