AM Recombinant Hyaluronidase

ARTSMedia Oocyte Denudation Medium

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ARTSMedia Denmark
Kongevejen 149,
DK-2830 Virum, Denmark

Indications for use

AM Recombinant Hyaluronidase is intended for denudation of oocytes prior to ICSI.

Device Description

1 ml in plastic bottles, individually packed, provided sterile. AM Recombinant Hyaluronidase is a HEPES buffered balanced salt solution, with 80 IU/ml recombinant hyaluronidase, bicarbonate, sodium pyruvate, minerals, synthetic serum replacement, human serum albumin (HSA) and gentamicin sulphate 10 µg/ml.

For professional use only.

Quality Control Specifications

pH: 7.2–7.6
(Ph Eur. 2.2.3, USP <791>)
Osmolality: 270-300 mOsm/kg
(Ph Eur. 2.2.35, USP <785>)

Appearance: Clear and particulate free
Sterility: No growth
(Ph Eur. 2.6.1, USP <71>)

Endotoxins: < 0.50 EU/ml (USP <85>)
1-cell Mouse Embryo Assay (MEA): ≥ 80% of embryos developed to expanded blastocysts at 96 hours (following 30 minutes exposure to AM Recombinant Hyaluronidase).
Enzymatic Activity: 40-160 IU/ml by the current USP test.

A certificate of analysis is available at www.artsmediadenmark.com

Storage

Store in the original container between 2-8 °C and avoid exposure to sunlight. Do not freeze.

Warnings

- Do not use after expiry date.
- Do not use if the product is cloudy or shows any evidence of microbial contamination.
- Do not use if the seal of the container is damaged or open upon delivery.
- Do not use medium that shows any evidence of particulate matter.
- Do not re-sterilize after opening.
- Do not re-use.
- Single use only.
- Discard excess (unused) media that has been warmed in accordance with local regulations for disposal of medi-

cal devices.

- AM Recombinant Hyaluronidase contains gentamicin sulphate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic (allergic reaction).

Precautions

- **Caution:** U.S. Federal law restricts this device to sale by or on the order of a physician.
- Media use should be conducted using aseptic techniques.
- Always wear protective clothing when handling samples.
- This product contains albumin, a derivative of human blood. All blood products should be treated as potentially infectious. Source material used for this product was tested and found negative for HBsAg, Anti-HIV-1/-HIV2, HCV, HCV-RNA, and HIV-RNA. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated.

General Notes Before Use

AM Recombinant Hyaluronidase is HEPES buffered and does not require CO₂ equilibration.

Instructions for Use

AM Recombinant Hyaluronidase should be used according to the established procedures and protocols optimized for the individual clinic.

Example

Preheat the AM Recombinant Hyaluronidase solution to 37 °C.

Prepare a drop of the preheated AM Recombinant Hyaluronidase solution on a culture dish or petri dish.

Transfer the cumulus oocyte complex (COCs) carefully into the drops for no longer one minute while using a denudation pipette.

Transfer the oocytes to a fresh dish with AM Handling Medium or AM-IVC Medium to wash away any excess hyaluronidase and dissociated cells. Gently pipette the oocytes in this fresh medium several times to ensure thorough removal of all residues, should any cumulus cells remain.

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