

## Material Safety Data Sheet (MSDS)

### ARTSMedia High Viscosity Oil (AM-HV Oil)

Product Number: 506050 and 506100 | Issue: 2 | Date of Revision: 2023DEC05

#### SECTION 1: Identification of the substance/mixture and of the company/undertaking

##### 1.1 Product identifier:

**Product form:** Mixture

**Product name:** ARTSMedia High Viscosity Oil (AM-HV Oil)

##### 1.2. Relevant identified uses of the substance or mixture and uses advised against

###### 1.2.1 relevant identified uses

**Industrial/professional use spec:** For professional use only

**Use of substance/mixture:** AM-HV Oil is a high viscosity mineral oil (d=0.86-0.88g/ml) pre-washed with ultra-pure water for use in In Vitro Fertilization (IVF), Intracytoplasmic Sperm Injection (ICSI) and assisted reproductive technologies (ART).

###### 1.2.2. Uses advised against

No additional information available

##### 1.3. Details of the supplier of the safety data sheet

ARTSMedia Denmark ApS

Kongevejen 149,

2830 Virum

Denmark

info@artsmedia.dk

**Customer service:** +45 5350 4530

##### 1.4. Emergency telephone number

**Emergency number:** +45 82121212

#### SECTION 2: Hazards identifications

##### 2.1. Classification of the substance or mixture

**Classification according to Regulation (EC) No. 1272/2008 [CLP]**

Not classified

**Adverse physicochemical, human health and environmental effects**

No additional information available

##### 2.2. Label elements

**Labelling according to Regulation (EC) No. 1272/2008 [CLP]**

No labelling applicable

##### 2.3. Other hazards

Contains no PBT/vPvB substances  $\geq 0.1\%$  in accordance with REACH Annex XIII.

The mixture does not contain substances included in the list established in accordance with article 59(1) of REACH for having endocrine disrupting properties or is not identified as having endocrine disrupting properties in accordance with the criteria set out in Commission Delegated Regulation (EU) 2017/2100 or commission Regulation (EU) 2018/605.

##### Company address:

ARTSMedia Denmark ApS  
Kongevejen 149  
2830 Virum  
Denmark

##### Company details:

CRN: 41808845  
VAT number: DK41808845  
Phone: +45 53 50 45 30  
info@artsmedia.dk

##### Bank account:

Bank: Sparekassen Sjælland  
Registration number: 6070  
Account number: 0000379518  
SWIFT: SWESDK22  
IBAN: DK3560700000379518

### SECTION 3: Composition/information on ingredients

#### 3.1. Substances

Not applicable

#### 3.2. Mixtures

Name	Product identifier	%	Classification according to Regulation (EC) No. 1272/2008 [CLP]
Mineral Oil	(CAS No.) 8042-47-5 (EC No.) 232-455-8	100	Not classified

Comments: ARTSMedia High Viscosity Oil does not contain any raw materials of direct animal-origin or materials that have been produced using animal-origin components. Components have not been in contact with material of animal-origin processing and therefore pose no TSE risk.

### SECTION 4: First aid measures

#### 4.1. Description of first aid measures

**First aid measures in general:** Never give anything by mouth to an unconscious person. If you feel unwell, seek medical advice (show the label where possible).

**First aid measures after inhalation:** Allow affected person to breath fresh air. Allow the victim to rest.

**First aid measures after skin contact:** Remove affected clothing and wash all exposed skin areas with mild soap and water, followed by warm water rinse.

**First aid measures after eye contact:** Rinse immediately with plenty of water. Obtain medical attention if pain, blinking or redness persists.

**First aid measures after ingestion:** Rinse mouth. Do NOT induce vomiting. Obtain emergency medical attention.

#### 4.2. Most important symptoms and effects, both acute and delayed.

Symptoms/effects: Not expected to present a significant hazard under anticipated conditions of normal use.

#### 4.3. Indication of any immediate medical attention and special treatment needed.

No additional information available.

### SECTION 5: Firefighting measures

#### 5.1. Extinguishing media

Suitable extinguishing media: Foam, Dry powder, Carbon dioxide, Water spray. Sand.

Unsuitable extinguishing media: Do not use a heavy water stream.

#### 5.2. Special hazards arising from the substance or mixture.

No additional information available

#### 5.3. Advice for firefighters

Firefighting instructions: Use water spray or fog for cooling exposed containers. Exercise caution when fighting any chemical fire. Prevent firefighting water from entering the environment.

Protection during firefighting: Do not enter the fire area without proper protective equipment, including respiratory protection.

#### Company address:

ARTSMedia Denmark ApS  
Kongevejen 149  
2830 Virum  
Denmark

#### Company details:

CRN: 41808845  
VAT number: DK41808845  
Phone: +45 53 50 45 30  
info@artsmedia.dk

#### Bank account:

Bank: Sparekassen Sjælland  
Registration number: 6070  
Account number: 0000379518  
SWIFT: SWESDK22  
IBAN: DK356070000379518

## SECTION 6: Accidental release measures

### 6.1. Personal precautions, protective equipment and emergency procedures

#### 6.1.1. For non-emergency personnel

Emergency procedures: Evacuate unnecessary personnel.

#### 6.1.2. For emergency responders

Protective equipment: Equip cleanup crew with proper protection.

Emergency procedures: Ventilate area.

### 6.2. Environmental precautions

Prevent entry to sewers and public waters.

### 6.3. Methods and material for containment and cleaning up

Methods for cleaning up: Soak up spills as soon as possible. Collect spillage. Store away from other materials.

### 6.4. Reference to other sections

See Heading 8. Exposure controls and personal protection.

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

Precautions for safe handling: Wash hands and other exposed areas with mild soap and water before eating, drinking or smoking and when leaving work. Provide good ventilation in process area to prevent formation of vapor.

### 7.2. Conditions for safe storage, including any incompatibilities.

Storage conditions: Keep only in the original container in a cool, well-ventilated place away from direct (sun)light. Keep container closed when not in use. Do not freeze. Do not use it after expiry date. After opening the container, the product can be used up to 1 month when sterile conditions are maintained, and the products are stored at room temperature.

Cannot be resterilized after opening. Do not expose product to sources of irradiation.

Incompatible products: Strong bases, strong acids.

Incompatible materials: Sources of ignition. Direct sunlight.

Storage temperature: 15-25°C

### 7.3. Specific end use(s)

See instructions for use delivered with the product.

## SECTION 8: Exposure controls/Personal protection

### 8.1. Control parameters

No additional information available

### 8.2. Exposure controls

#### Appropriate engineering control:

Handle in accordance with good industrial hygiene and safety. Avoid all unnecessary exposure.

#### Personal protective equipment:

Wear fire/flame resistant/retardant clothing.

#### Skin protection:

Complete suit protecting against chemicals, flame retardant antistatic protective clothing. The type of protection must be selected according to the concentration and volume of the dangerous substance at the specific workplace.

#### Hand protection:

#### Company address:

ARTSMedia Denmark ApS  
Kongevejen 149  
2830 Virum  
Denmark

#### Company details:

CRN: 41808845  
VAT number: DK41808845  
Phone: +45 53 50 45 30  
info@artsmmedia.dk

#### Bank account:

Bank: Sparekassen Sjælland  
Registration number: 6070  
Account number: 0000379518  
SWIFT: SWESDK22  
IBAN: DK3560700000379518

Wear protective gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.

The selected protective gloves have to meet specifications of regulation EU 2016/425 and the standard EN374.

**Eye/Face protection:**

Face shield and safety glasses. Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166 (EU).

**Respiratory protection:**

Wear an appropriate mask. Where risk assessment shows air-purifying respirators are appropriate, use a full-face respirator with multi-purpose combination (US) or type ABEK (EN 14387) respirator cartridges as a backup to engineering controls. If the respirator is the sole means of protection, use a full-face supplied air respirator. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

**Other information:**

Do not eat, drink or smoke during use. Do not pipette liquid using a mouth pipette.

**SECTION 9: Physical and chemical properties****9.1. Information on basic physical and chemical properties**

Physical state	:	Liquid
Colour	:	Transparant liquid
Apperance	:	Clear, colorless liquid
Odour	:	Odourless
Odour threshold	:	No data available
pH	:	No data available
Melting point	:	No data available
Freezing point	:	No data available
Boiling point	:	No data available
Flash point	:	No data available
Auto-ignition temperature	:	No data available
Decomposition temperature	:	No data available
Flammability	:	Non flammable.
Vapour pressure	:	No data available
Vapour pressure at 50°C	:	No data available
Relative vapour density at 20°C	:	No data available
Density	:	No data available
Relative density	:	0.86-0.88 g/cm <sup>3</sup> @ 25°C
Solubility	:	Not soluble in water
Partition coefficient n-octanol/water (Log Kow)	:	No data available
Viscosity, kinematic	:	No data available
Particle characteristics	:	Not applicable
Explosive limits	:	No data available
Viscosity	:	95-140 cP @ 30°C

**9.2. Other information**

No additional information available

**SECTION 10: Stability and Reactivity****Company address:**

ARTSMedia Denmark ApS  
Kongevejen 149  
2830 Virum  
Denmark

**Company details:**

CRN: 41808845  
VAT number: DK41808845  
Phone: +45 53 50 45 30  
info@artsmedia.dk

**Bank account:**

Bank: Sparekassen Sjælland  
Registration number: 6070  
Account number: 0000379518  
SWIFT: SWESDK22  
IBAN: DK3560700000379518

**10.1. Reactivity**

No additional information available

**10.2. Chemical stability**

Stable after transport (max. 5 days) at elevated temperature ( $\leq 37^{\circ}\text{C}$ ). Aim stability for 18 months from date of manufacture.

**10.3. Possibility of hazardous reactions**

Not established.

**10.4. Conditions to avoid.**

Direct (sun)light. Extremely high or low temperatures. Do not expose products to sources of irradiation.

**10.5. Incompatible materials**

Strong acids. Strong bases.

**10.6. Hazardous decomposition products**

Carbon monoxide. Carbon dioxide. fume.

**SECTION 11: Toxicological information****11.1 Information on toxicological effects**

**Acute toxicity (oral /dermal/inhalation):** Not classified.

Extensive data on Mouse Embryo Assays have demonstrated that the medium is not toxic.

**Skin corrosion/irritation:** Not classified. Based on available data, the classification criteria are not met

**Serious eye damage/irritation:** Based on available data, the classification criteria are not met.

**Respiratory or skin sensitisation:** Not classified. Based on available data, the classification criteria are not met

**Germ cell mutagenicity:** Not classified. Based on available data, the classification criteria are not met.

**Carcinogenicity:** Not classified. Based on available data, the classification criteria are not met

**Reproductive toxicity:** Not classified. Based on available data, the classification criteria are not met

**STOT-single exposure:** Not classified. Based on available data, the classification criteria are not met.

**STOT-repeated exposure:** Not classified. Based on available data, the classification criteria are not met.

**Aspiration hazard:** Not classified. Based on available data, the classification criteria are not met.

**11.2 Information on other hazards:**

**Endocrine disrupting Properties:** Based on available data, the classification criteria are not met

**Potential adverse human health effects and symptoms:** Based on available data, the classification criteria are not met.

**Additional information:** Based on available data, the classification criteria are not met.

**SECTION 12: Ecological information****12.1. Toxicity**

Hazardous to the aquatic environment, short-term (acute): Not classified

Hazardous to the aquatic environment, long-term (chronic): Not classified

**12.2 Persistence and degradability**

Persistence and degradability

Not established.

**Company address:**

ARTSMedia Denmark ApS  
Kongevejen 149  
2830 Virum  
Denmark

**Company details:**

CRN: 41808845  
VAT number: DK41808845  
Phone: +45 53 50 45 30  
info@artsmedia.dk

**Bank account:**

Bank: Sparekassen Sjælland  
Registration number: 6070  
Account number: 0000379518  
SWIFT: SWESDK22  
IBAN: DK3560700000379518

**12.3 Bioaccumulative potential**

Bioaccumulative potential

Not established.

**12.4 Mobility in soil**

No additional information available

**12.5 Results of PBT and vPvB assessment**

No additional information available

**12.6 Endocrine disrupting Properties**

No additional information available

**12.6. Other adverse effects**

Avoid release to the environment.

**SECTION 13: Disposal considerations****13.1. Waste treatment methods****Product/Packaging disposal recommendations:**

Dispose in a safe manner in accordance with local/national regulations.

**Ecology - waste materials:**

Avoid release to the environment.

**SECTION 14: Transport information**

In accordance with ADR / RID / IMDG / IATA / ADN

**14.1. UN number**

UN-No. (ADR/ IMDG/IATA/ADN/RID): Not applicable

**14.2. UN proper shipping name****Proper Shipping Name**

(ADR/ IMDG/IATA/ ADN /RID): Not applicable

**14.3. Transport hazard class(es)****Transport hazard class(es)**

(ADR/ IMDG/IATA/ ADN /RID): Not applicable

**14.4. Packing group**

Packing group (ADR/ IMDG/IATA/ ADN /RID): Not applicable

**14.5. Environmental hazards****Dangerous for the environment:** No**Marine pollutant:** No**Other information:** No supplementary information available**14.6. Special precautions for user****Overland transport:** No data available**Transport by sea:** No data available**Air transport:** No data available**Inland waterway transport:** No data available**Rail transport:** No data available**14.7. Transport in bulk according to Annex II of Marpol and the IBC Code:** Not applicable**Company address:**ARTSMedia Denmark ApS  
Kongevejen 149  
2830 Virum  
Denmark**Company details:**CRN: 41808845  
VAT number: DK41808845  
Phone: +45 53 50 45 30  
info@artsmmedia.dk**Bank account:**Bank: Sparekassen Sjælland  
Registration number: 6070  
Account number: 0000379518  
SWIFT: SWESDK22  
IBAN: DK3560700000379518

## SECTION 15: Regulatory information

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

#### 15.1.1. EU-Regulations

Contains no REACH substances with Annex XVII restrictions

Contains no substance on the REACH candidate list

Contains no REACH Annex XIV substances

Contains no substance subject to REGULATION (EU) No 649/2012 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 4 July 2012 concerning the export and import of hazardous chemicals.

Contains no substance subject to Regulation (EU) No 2019/1021 of the European Parliament and of the Council of 20 June 2019 on persistent organic pollutants

### 15.2. Chemical safety assessment

No chemical safety assessment has been carried out

## SECTION 16: Other information

### Data sources:

REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

**Other information:** None.

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product.

#### Company address:

ARTSMedia Denmark ApS  
Kongevejen 149  
2830 Virum  
Denmark

#### Company details:

CRN: 41808845  
VAT number: DK41808845  
Phone: +45 53 50 45 30  
info@artsmedia.dk

#### Bank account:

Bank: Sparekassen Sjælland  
Registration number: 6070  
Account number: 0000379518  
SWIFT: SWESDK22  
IBAN: DK3560700000379518